

Photo Rejuvenation

Hair Removal

Photo Facial

Acne

Vascular Lesions

Large Pores


Rosacea

Freckles

Brown Spots

Sun Damaged Skin

Skin Tightening



5025 Orbitor Drive
Building 6, Suite 301
Mississauga, Ontario
L4W 4Y5, Canada
Tel: 416.628.7191
Fax: 416.628.9940
www.cmedtech.ca




SOLISPRO SYSTEM™
Photo Rejuvenation & Hair Removal

Photo Facial



Hair Removal




Intelligent Contact Cooling completely protects the epidermis up to 0.5mm below the skin. With a maximum cooling temperature down to -10 C for a safe, effective and painless treatment.

Photo Rejuvenation



Acne



Integrated Temperature Sensor that reads the temperature of the skin during treatment and only allows emission when the specified skin temperature is reached.

Vascular Lesions



Sun Damaged Skin

