


Acne Scars
Skin Tightening
Wrinkle Reduction
Pore Reduction
Skin Rejuvenation
Scars
Non-Surgical Facelift



No Downtime
No Pigment Change



Jeisys

The Leader in Medical Technology

5025 Orbitor Drive, Building 6, Suite 301
Mississauga, Ontario, L4W 4Y5, Canada
Tel: 416.628.7191

WWW.CMEDTECH.CA

ACNE SCARS



Before Treatment



34 Days After Treatment

WRINKLE REDUCTION




Before Treatment




58 Days After Treatment

HOW INTRAcel WORKS



Tissue Damage - Fractional Laser vs INTRAcel



The INTRAcel RF Microneedles are able to apply the energy directly to the target area. Fractional Lasers require high energy levels to attempt to reach the target area causing significant surface damage to the skin as well as pigment changes and creating longer recovery time after treatment.